



Infrared Controlled Telephones

An infrared controlled telephone can be used by virtually any Electronic Daily Living Aid (EADL, formally Environmental Controls) that can send infrared signal. Here is a comparison of features:

Features	GewaTel 200 ZYGO Industries, Inc. 800-234-6006 www.zygo-usa.com	Infrared Home or Office Accessible Landline Telephone Broadened Horizons 612-851-1040 www.broadenedhorizons.com
Photo		
cost	\$1,540.00	\$699.00
access	-IR signal -Keypad	-IR signal -Keypad -Voice with voice dialer
user can build #s	Yes	Yes
prestored #s held	50	50
redial	Yes	Yes
call waiting	Yes	Yes
battery back-up	Yes	Yes
Display	Yes	Yes
Volume Control	Yes	Yes
Privacy headset	\$115 earphone \$340 microphone	Wired privacy headset \$29.95 Wireless headset \$349.95 Expansion Microphone \$149.95

Interfaces	-any learning IR transmitter -any IR transmitting EADL -GEWA Prog	-any learning IR transmitter -any IR transmitting EADL -VoiceIR Voice Controller
Transmission	Analog landline	Analog landline
Other	-can be used as standard phone.	-can be used as a standard phone -voice dialer \$249.95 -same base phone as GewaTel 200

The GEWATel200 and Infrared Home or Office Accessible Landline Telephone are a Konftel 200 phone with remote. These are often available through various internet retailers for a lesser cost.

Michelle L. Lange, OTR, ABDA, ATP, Access to Independence, MichelleLange@msn.com, 7/10.